
lisa


Converter *lisa* → Funk
Converter *lisa* → radio
Convertiseur *lisa* → radio
Convertor *lisa* → radio
Convertitore *lisa* → radio
Convertidor *lisa* → radio

Bedienungsanleitung
Operating Instructions
Mode d'emploi
Gebruiksaanwijzing
Istruzioni per l'uso
Instrucciones de servicio

Vorderseite

Rückseite

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

Front

Back

Congratulations on purchasing your *lisa* → *Radio Converter*. We hope that it will serve you well for many years to come. You have opted for a modern and reliable system. Please read through these operating instructions carefully. They tell you how to set up the system correctly and explain all the system's features to you.

Functional principle

With the *lisa* → *Radio Converter*, *lisa* systems operating through the 230 V power supply can be combined with the radio receivers. The *lisa* → *Radio Converter* receives the signals from various *lisa* transmitters and sends them to the cable-free radio receivers. All *lisa* transmitters (acoustic/galvanic telephone transmitters, acoustic/galvanic door bell transmitters, baby phones, alarm transmitters) can be used as transmitters.

Operating and display elements

- Person call button
- Display LEDs
- Channel selection switch
- Connection socket for vibration cushion
- Transmit with delay (this operating element is only important if two *lisa* radio converters are used).

Note 1


The person call button located on the converter does not have any function for 60 seconds if an alarm signal/person call signal has been sent.

Note 2

If the "Transmit with delay" switch is on position 1 (pressed in) the radio signals are sent with a delay of 1.5 seconds.

Relay function – operation with two *lisa* → *Radio Converters*

If the radio transmission range of a *lisa* → *Radio Converter* is not sufficient (e.g. in large houses), you can incorporate a second *lisa* → *Radio Converter* into the system.

Important: *lisa* → radio converter 1 :
 "Transmit with delay" is pressed


lisa → radio converter 2 :
 : "Transmit with delay" is not pressed

Channel selection

A digital transmission system ensures the correct recognition of each of the digital signals from each transmitter. The system has 10 different selectable channels, 0-9. The factory setting for all *lisa* devices is channel 0. Only where more than one installation is required in the same building will it be necessary to select a different channel for each system.

To select a channel, turn the rotary switch on the underside of the transmitter to the required channel number using a small screwdriver.

Example:
 Set channel number 0

Startup

The *lisa* → *Radio Converter* is immediately ready for use after being connected to a 230V socket.

Transmission range

The transmission range of the radio signals of the *lisa* → *Radio Converter* is designed for 80 metres under optimal conditions.

The following conditions may result in a reduced range:

- Signal transmission through building structures or vegetation
- Interference from TVs, computer monitors and mobile phones. Under certain circumstances, this may prevent the receiver from functioning at all.
- Positioning of the transmitter at an unfavorable distance from surfaces which reflect radio waves, for example floors and walls. This may lead to the signal wave being weakened or even cancelled out by the reflected wave.
- Metallic objects reduce the range due to the shielding effect of the metal.

- Particularly in towns and cities, there may be many other sources of radio waves which can disrupt the original signal.
- Having units operating at similar frequencies located close to one another may also result in mutual interference.

Positioning

In order to attain as wide a range as possible, you should follow the installation information below :

- The *lisa* → *Radio Converter* should not be positioned near devices which emit strong interference signals such as television sets, monitors or mobile telephones.
- The *lisa* → *Radio Converter* should not be positioned on metallic objects as otherwise the transmission range is reduced by the metallic screening.
- Position the *lisa* → *Radio Converter* as centrally as possible in the house and at a medium room height.

Maintenance and care

All *lisa* devices are maintenance free. To clean, remove from their mains sockets and wipe over with a moistened cloth. Never use alcohol, thinner or other organic materials.

lisa devices should not be exposed to excessive direct sunlight, heat, thick dust, humidity or mechanical shock.

Warranty

This warranty covers the repair of the product and returning it to you free of charge. It is essential that you send in the product in its original packaging, so do not throw the packaging away.

The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people who are not authorised to do so (destruction of the seal on the unit). Repairs under warranty are only carried out providing you have filled in and returned the enclosed warranty card from the dealer.

Technical data

Power supply: 230 V, 50 Hz
Frequency: 868.35 MHz
Power consumpt. approx. 1 W


Height: 40 mm
Width: 120 mm
Depth: 135 mm
Weight: approx. 350 g

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Humantechnik GHL AG
ELECTRONIQUE S.A.R.L