

lisa ((•))

D	Bedienungsanleitung Funk-Converter	Seite 2
GB	Operating Instructions Radio-converter	Page 6
F	Mode d'emploi Convertisseur radio	Page 10
NL	Gebruiksaanwijzing Converter radio signaal	Pagina 14
I	Istruzioni per l'uso Convertitore radio	Pagina 18
E	Instrucciones de servicio Convertidor radio	Página 22

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
Congratulations on purchasing your *Radio-lisa converter*. We hope it will give you many years of reassuring operation. You have purchased a modern and reliable system. Please read through these operating instructions carefully. They tell you how to set up the system correctly and explain all the system's features to you.

Functional principle

With the *Radio-lisa converter* it is possible to combine the *lisa* systems, which is working on the 230V-powerline, with the Radio Transmitters. The *Radio-lisa converter* receives the signals from different Radio Transmitters and feed this signals into the 230V-powerline. Following radio transmitters can be used: *Radio Telephone Transmitter acoustic/ direct connection*, *Radio Doorbell Transmitter acoustic/ direct connection*, *Radio Baby Monitor*, *Radio Alarm Transmitter*.

Controls


- » Transmitter« indicator
- Call button
- Channel selector

Note: The pager button on the *Radio-lisa converter* ceases to function for 40 seconds once a manual or alarm call has been detected.

Channel selection

A digital transmission system ensures the correct recognition of each of the digital signals from each transmitter. The system has 10 different selectable channels, 0-9. The factory setting for all *lisa* devices is channel 0. Only where more than one installation is required in the same building will it be necessary to select a different channel for each system.

To select a channel, turn the rotary switch on the underside of the transmitter to the required channel number using a small screwdriver.

All transmitters and receivers must be set to the same channel.

Installing the *Radio-lisa converter*

Simply connect the *Radio-lisa converter* to a power point. The unit is immediately ready to use.

Radio range

The transmitter for the *Radio-lisa converter* are designed to have a radio signal range of 80 meters under optimum conditions.

The following conditions may result in a reduced range:

- Signal transmission through building structures or vegetation
- Interference from TVs, computer monitors and mobile phones. Under certain circumstances, this may prevent the receiver from functioning at all.
- Positioning of the transmitter at an unfavourable distance from surfaces which reflect radio waves, for example floors and walls. This may lead to the signal wave being weakened or even cancelled out by the reflected wave.

- Metallic objects reduce the range due to the shielding effect of the metal.
- Particularly in towns and cities, there may be many other sources of radio waves which can disrupt the original signal.
- Having units operating at similar frequencies located close to one another may also result in mutual interference.

Positioning the radio transmitter

Bear in mind the following points when positioning the radio transmitter, in order to achieve the greatest possible range from the radio transmitter:

- The radio transmitter should not be set up close to devices which emit powerful interference signals, such as TVs, computer monitors or mobile phones.
- The radio transmitter should not be placed on metallic objects, since otherwise the radio range will be restricted by metallic shielding.
- Position the radio transmitter as close to the centre of the building as possible, at a medium height in the room.

Maintenance and care

All *lisa* devices are maintenance free. To clean, remove from their mains sockets and wipe over with a moistened cloth. Never use alcohol, thinner or other organic materials.

lisa devices should not be exposed to excessive direct sunlight, heat, thick dust, humidity or mechanical shock.

Warranty

This warranty covers the repair of the product and returning it to you free of charge. It is essential that you send in the product in its original packaging, so do not throw the packaging away.

The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people who are not authorised to do so (destruction of the seal on the unit).



Repairs under warranty are only carried out providing you have filled in and returned the enclosed warranty card from the dealer.

Technical data

Power supply: 230V/ 50Hz
Frequency: 868.35 MHz
Operating range: 0 – 40 °C

Height: 21 mm
Width: 65 mm
Depth: 35 mm
Weight: approx. 200 g

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Humantechnik GHL AG