

<i>lisa / lisa RF</i>	Sender
<i>ISDN /</i>	Transmitter
<i>Numéris</i>	Émetteur
	Zender
	Trasmittitore
	Emisor

Bedienungsanleitung
Operating Instructions
Mode d'emploi
Gebruiksaanwijzing
Istruzioni per l'uso
Instrucciones de servicio

Allgemeines

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

General

The *lisa* ISDN transmitter / *lisa* RF ISDN transmitter has been specially developed for connection to the Euro-ISDN network.

It reacts to an incoming phone call for one or more ISDN terminal units. When a phone call is detected, radio pulses are transmitted to the receivers via the standard socket (triggering a phone signal). The receiver converts the radio pulses into light signals/vibrations.

The three different LEDs on the *lisa* ISDN transmitter / *lisa* RF ISDN transmitter enable you to tell which ISDN terminal unit is ringing.

Note

The transmitter does not react if no terminal unit is connected to your ISDN phone connection. The same applies if the ringer of your terminal unit is deactivated.


Setting the correct radio channel

A special digital transmission process ensures that the radio signals are picked up without errors. This process also makes it possible to set 10 different radio channels. When supplied from the factory, all *lisa* units are set to channel 0.

It is only necessary to set units to other radio channels if there are several systems being operated at the same time in the same building.

How to set the channel

Use a small screwdriver to turn the channel selector switch on the back of the unit until the arrow is pointing to the required channel number.

Example:
Set channel
number 0

Pager button

Pressing the pager button triggers a certain sequence of flashes/vibrations on every receiver. This allows you to call people in other rooms. Also, this button can be used for checking the system is functioning even when no corresponding cables are connected.

**Changing the flash duration
(Lisa ISDN transmitter only!)**

When new, all units are set so each Lisa receiver flashes for 28 seconds. However, you can reduce this time to 8 seconds.

Procedure

Disconnect the mains lead of the transmitter from the power socket and reconnect it whilst holding the pager button pressed. This sets the flash duration to 8 seconds. If you want to restore the flash duration to 28 seconds, simply disconnect and reconnect the mains lead of the transmitter to the wall socket (this time without pressing the pager button).

Connecting the ISDN transmitter

The ISDN transmitter has two ISDN sockets. These sockets are connected in parallel inside the unit. This makes it very simple to connect to ISDN terminal units.

Procedure

Remove the ISDN connection cable from your ISDN terminal unit (e.g. telephone) and connect it to the back of the ISDN transmitter. Using the supplied ISDN cable, connect the ISDN transmitter to your ISDN terminal unit (e.g. telephone).

Delivery condition

No phone numbers (MSNs) are stored in the ISDN transmitter when it is new. The "P" LED flashes to indicate that the unit is not programmed.

Unprogrammed status

The LED marked "P" flashes. Every incoming phone call, irrespective of which phone number (MSN) is dialled, triggers a transmit signal. The LED marked "1" is displayed on the ISDN transmitter for each phone call.

Programming the ISDN transmitter

The ISDN transmitter incorporates innovative technology which makes programming much easier. Simply use another phone or mobile phone to call the ISDN terminal unit which you want to monitor with the ISDN transmitter.

The ISDN transmitter has three memory slots which can be used for monitoring three different phone numbers (MSNs).

Procedure

- 1 Press the button marked "PR" (hold it pressed for one second).
- 2 The "P" LED flashes steadily to indicate the unit is in programming mode.
- 3 You can use the "PR" button to select one of the three memory slots of the ISDN transmitter. Only press the "PR" button briefly to do this. The current position is shown by a flashing light.
Note: Already programmed memory slots are illuminated. It is also possible to overwrite a programmed memory slot.
- 4 Dial the ISDN terminal unit which you want to monitor. When the ISDN terminal unit rings, the phone number (MSN) is automatically transferred into the ISDN transmitter. The ISDN transmitter exits programming mode if the programming was successful. Repeat the procedure for each of ISDN terminal unit which you want to monitor (maximum 3).

Important

- As soon as one phone number has been programmed in the ISDN transmitter, the transmitter only responds to programmed phone numbers (MSNs).

- Phone numbers which have already been stored can only be overwritten. Selective deletion of individual memory slots is not possible.

Reset

Press the »R« button on the back of the unit for 5 seconds. All programmed phone numbers (MSNs) are deleted. The »P« LED on the front flashes (unprogrammed status).

Maintenance and care

The *lisa* ISDN transmitter / *lisa* RF ISDN transmitter is maintenance-free. If the unit does become dirty, simply wipe it clean with a soft, damp cloth. Never use spirits, thinners or other organic solvents. Do not set up the unit where it will be exposed to full sunlight for long periods. In addition, it must be protected against excessive heat, moisture and severe mechanical shocks.

Warranty



The *lisa* ISDN transmitter / *lisa* RF ISDN transmitter is highly reliable. However, please contact your dealer or the manufacturer directly if a malfunction occurs in spite of the unit having been set up and operated correctly. This warranty covers the repair of the product and returning it to you free of charge. It is essential to return the phone in its original packaging, so do not throw this away. The warranty does not cover damage caused by incorrect handling or attempted repairs by unauthorised persons (destruction of the seal on the unit). Repairs will only be carried out under warranty if the dealer's warranty card has been filled in and returned or if you have the original warranty certificate if the phone was purchased directly from the manufacturer. The serial number must be specified in all cases.

Technical data

Power supply: 230 V, 50 Hz
Power consumpt. approx: 1 W

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Humantechnik GHL AG
ELECTRONIQUE S.A.R.L