

*lisa*

|           |  |           |
|-----------|--|-----------|
| <b>D</b>  | <b>Bedienungsanleitung</b><br><i>lisa time S</i>       | Seite 2   |
| <b>GB</b> | <b>Operating Instructions</b><br><i>lisa time S</i>    | Page 6    |
| <b>F</b>  | <b>Mode d'emploi</b><br><i>lisa time S</i>             | Page 10   |
| <b>NL</b> | <b>Gebruiksaanwijzing</b><br><i>lisa time S</i>        | Pagina 14 |
| <b>I</b>  | <b>Istruzioni per l'uso</b><br><i>lisa time S</i>      | Pagina 18 |
| <b>E</b>  | <b>Instrucciones de servicio</b><br><i>lisa time S</i> | Página 22 |

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

***lisa time S*** is an electronic digital clock with a receiver to the *lisa* signalling system. It wakes you up using a flashing light and/or vibration and/or acoustic signal. Additionally, the wake-up or *lisa* signal is issued through a switched socket.


#### Using your clock

Plug the mains plug of the *lisa time S* into a 230V mains socket. The display flashes as a sign that the time is not correctly set. As soon as the time has been set the flashing stops.


#### Setting the time

Keep the button  and the **h** or **min** button pressed simultaneously to set the desired time.

#### Setting the alarm time

Keep the button  and the **h** or **min** button pressed simultaneously to set the desired alarm time.

#### Checking the alarm time set

Press the button , the alarm time shown is in the display of the alarm clock. A flashing dot between the hour

and the minutes on the display shows that the alarm has been set.

#### Alarm signal

After reaching the alarm time set the *lisa time S* starts to issue flashing lights and/or an acoustic signal depending on the setting. If a vibration cushion is connected the vibration pulse is also simultaneously triggered. With the switched socket the alarm signal can be passed to a lamp or a vibration system. If an alarm signal is not interrupted by pressing the "Alarm Stop 9 min" or the "Alarm Stop 24h" button, the alarm signal is repeated after 4 minutes (auto snooze).

#### Interrupting/switching off the alarm signal:

##### For 9 minutes (snooze):

By pressing the snooze button "Alarm Stop 9 min"

##### For 24 hours:

By pressing the "Alarm Stop 24 h"

##### For a longer period (week-end, holiday):

Reset the switch on the back to **0**

### **Special features of *lisa time S***

- When using the alarm clock for the first time it is possible that the internal battery needs recharging. This means that the alarm clock will only be fully functional after about six hours.
- *lisa time S* is also a receiver for the *lisa* signalling system. From the 4 signal lights you can recognise which alarm signals (door, telephone, baby, alarm or beeper) are being issued. As with all *lisa* devices *lisa time S* has a switch to change channels (also see *lisa* operating instructions).
- The brightness of the display can be set by a potentiometer.
- The built-in battery can bridge a power failure of at least one hour. During this time the time and alarm setting remain saved.
- If a light is connected it can be switched on when the switch is in the "light on" position. Please note that when the switch is in this position the switch output does not have an alarm function.

### **Warranty**

This warranty covers the repair of the product and returning it to you free of charge. It is essential that you send in the product in its original packaging, so do not throw the packaging away.



The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people who are not authorised to do so (destruction of the seal on the unit). Repairs under warranty are only carried out providing you have filled in and returned the enclosed warranty card from the dealer.

**Technical data**

|                                     |              |
|-------------------------------------|--------------|
| Power supply:                       | 230 V, 50 Hz |
| Power consumption:                  | approx. 1 W  |
| Maximum load of the switched socket | 500 W        |

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**Humantechnik GHL AG**